

**CERTIFICATE OF COMPLIANCE FOR RADIOACTIVE MATERIALS PACKAGES
D.O.T. SPECIFICATIONS 7A, TYPE A**

**Model #001-787 PET Shipping System
Page 1 of 2**

1. Package Identification:

Common Name: Double Pig Shipping System Syringe w/Needle
Ground or Air Shipment
Catalog Name: **Shipping System, PET, Double Pig, w/Needle**
Model Number: **#001-787**

2. Preamble:

- a. This certificate is issued to satisfy section 178.350 including sections 171.2, 173.403, 173.410, 173.412, 173.415, 173.431, 173.465, 173.466, 178.350 and 178.608 of the Department of Transportation Hazardous Materials Regulations (49 CFR).
- b. The packaging and contents described in item 4 below were found to meet the safety standards set forth in, and tested, according to the more stringent recommendations for general performance testing from:

IATA Dangerous Goods Regulations, 59th Edition Sections 5.0.4.3, 10.5, and 10.6.1 thru 10.6.3.5.

49 CFR as listed in item 2a above.
- c. This certificate does not release the consignor from compliance with any requirement of the regulations of the U.S. Department of Transportation or other applicable regulatory agencies, including the government of any country through or into which the package will be transported.

3. Descriptions of Packaging and Authorized contents, Other Conditions, and References:

- a. **Outer Container:** Model #001-783-A005
Medium Density Polyethylene:
Dimension: 11.75" x 11.75" x 12.5" with high density polyurethane foam cushioning and insert support weldment with two locking twist hasps and a hinged lid. Case has two handles. A safety seal can be put thru the hasp and the stand off for tamper proof evidence.
- b. **Shielded Insert:** Model #001-782-C200
Contour lead shield with open top and bottom and more lead in the vertical center. The lead inner shield thickness varies from 0.06 inches to 1.3 inches.
- c. **PET Unit Dose Pig:** Model #001-785 (System accepts 2 pigs)
0.065" polycarbonate with .5" lead with single turn twist lock top with 1.44" top lead and 1.2" lead in bottom, and an integral plastic liner or use pig Model #001-777 or #001-798 without liner – fits 10ml syringe with needle
- d. **Absorbent Material:** Model #001-771
One piece of 3" x 3" Quick Solid absorbent material in the bottom of each pig
- e. **Syringe:** B-D or Monoject plastic 3 ml or 5/6 ml Syringe with Needle

**CERTIFICATE OF COMPLIANCE FOR RADIOACTIVE MATERIALS PACKAGES
D.O.T. SPECIFICATIONS 7A, TYPE A**

**Model #001-787 PET Shipping System
Page 2 of 2**

Authorized Contents:

Normal Form Radioactive Materials in Type A quantities
Tested with F-18, 235 mCi and 160 mCi in appropriate parts of the shield as Yellow II.
Must comply with A₂ value limits from 49CFR 173.435.

Specifications and Restrictions:


Ground or Air Transportation

4. This certificate is issued on the basis of the test results and safety analysis of the package design.

5. Certification

This document is to certify that to the best of our knowledge, information and belief, the construction methods, the packaging design, and the materials of construction for this specification package comply with the standards described in items 2a, 2b, 2c herein. Any deviation from these conditions must be approved prior to use.

Approved:



Clyde Schlein
Vice President, Regulatory Affairs & Compliance
Biodex Medical Systems, Inc.

Date: Jan 4, 2013
Rev: Nov 11, 2015
Rev: Dec 5, 2017
Rev: Aug 14, 2018