

"The Clinical Advantage"™

Biodex Harnesses

For Biodex FreeStep SAS and Unweighing System



Unloading & Support

Patient comfort, safety and support for use on the Biodex FreeStep SAS or Biodex Unweighing System.

Harnesses are designed to accommodate a range of patient populations.

www.biodex.com/harnesses

NEW HARNESS

Biodex Standard Support Harness helps prevent slippage without limiting lower extremity range of motion.

BIODEX

www.biodex.com
1-800-224-6339

Int'l 631-924-9000

BIODEX HARNESSSES - UNLOADING & SUPPORT

FOR BIODEX FREESTEP SAS AND UNWEIGHING SYSTEM

NEW Standard Unweighing Harness

Designed for patient "unloading", provides safety and support
– This style is most popular for Biodex Unweighing System



The Biodex Unweighing Harness is designed to provide vertical lift through three shared points of pressure: through the hip straps, gluteal fold straps and the lumbar-thoracic harness. Its secure fit with maximum support offers comfort without limiting lower extremity range of motion. Assuring that the patient is comfortable is essential when lengthy exercise sessions are required.

The two-piece panel design and detachable pelvic support piece allow for a wide range of adjustability and easy application.

SPECIFICATIONS:

- Torso Circumference: adjustable from 24"-54" (61 to 137 cm)
- User Capacity: 300 lb (136 kg)
- Warranty: 60 days, parts

950-487 Harness, Standard Unweighing
For use with *Biodex Unweighing System*
and *FreeStep SAS*
NOTE: Standard Unweighing Harness included with Biodex Unweighing System

Watch Video: [YouTube www.biodex.com/harnesses-video](https://www.youtube.com/watch?v=www.biodex.com/harnesses-video)

Pediatric Harness

Support for the lightweights



For use with the Biodex Unweighing System or FreeStep SAS, the Pediatric Harness supports and unloads for children-size torsos.

SPECIFICATIONS:

- Torso Circumference: 14" to 36" (36 to 91 cm)
- Patient Capacity: rated for 30 to 100 lb (13 to 45 kg)
- Warranty: 60 days, parts

960-402 Harness, Pediatric
For use with *Biodex Unweighing System*
and *FreeStep SAS*

Support Harness DLX

A comfortable, supportive harness that is quick and easy to put on and off.



Worn as a vest around the upper torso, the Support Harness DLX provides patients with security and safety from falling when used with the FreeStep SAS or Unweighing System. Versatility includes double shoulder point connections to spreader bar or, a single connection point for low ceiling applications with tall patients (6'3" on a treadmill or 6'11" on the floor). A pelvic support strap is included to keep the harness in place when limited body weight support is intended.

Note: This vest is not intended for patient "unloading", but can be used with the Biodex Unweighing System when support is required.

SPECIFICATIONS:

- Torso Circumference: 28" to 50" (71 to 127 cm)
- Patient Capacity: rated for 300 lb (136 kg)
- Warranty: 60 days, parts

960-126 Harness, Support, DLX
Includes *Pelvic Support Strap*
For use with the *Biodex Unweighing System*
and *FreeStep SAS*

Bariatric Harness

Support for larger patients



Similar in design to the Standard Unweighing Harness, the Bariatric Harness grips patients around the torso and thighs for a secure fit when offloading a patient's weight during partial weight-bearing therapy (PWBT) on the Biodex Unweighing System.

SPECIFICATIONS:

- Torso Circumference: 52" to 72" (132 to 183 cm)
- Patient Capacity: rated for 600 lb (272 kg)
- Warranty: 60 days, parts

960-401 Harness, Bariatric
For use with *Biodex Unweighing System*
and *FreeStep SAS*

www.biodex.com/harnesses

BIODEX

Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9338, Email: info@biodex.com, www.biodex.com