


# UPPER BODY CYCLE



*EZ Grip Handgrips*  
accommodate patients with  
hand control limitations

*Tilting Display* for best viewing angle  
whether seated or standing

Easy-to-use "*Quick Start*" control  
panel - with large buttons

*Removable Seat* for  
wheelchair or  
standing access

*Adjustable Arm Cranks*  
to limit/increase range  
of motion

*Pivoting Actuator*  
accommodates a wide  
range of cycling positions

*Large Rotating Seat* for  
patients with  
ambulatory difficulties

*Grab Handles* provide  
added security

*Seat Belt* for torso  
stabilization



*Seat Front to Back Positioning Scale*  
for reference point

*Large Wheels* make  
transporting &  
repositioning easy

[www.biodex.com/ubc](http://www.biodex.com/ubc)

## BIODEX

Biodex Medical Systems, Inc.

20 Ramsey Road, Shirley, New York, 11967-4704, Tel: 800-224-6339 (Int'l 631-924-9000), Fax: 631-924-9241, Email: info@biodex.com, www.biodex.com